



**High acceleration.
Low jitter.
FASTER TESTING.**

The all-new DMM-101 Single Axis Controller delivers unmatched velocity stability performance, providing a stable test platform for increased aerial track densities. The DMM 101's flexible connector interface accepts both MCS-LA2000 and MFM-BDC-610 connections allowing it to be a drop-in replacement in today's media testers.

Our all-new ST1 Air Bearing Spindle offers extremely low asynchronous error motion (AEM) performance. Fast acceleration of less than 1 second to 15,000 RPM for 65 mm disks delivers improved throughput and faster testing of media.

Combine our DMM-101 controller and our ST1 spindle for a matched, easily integrated, sub-system solution for HDD applications.

Contact Danaher Motion - Dover for additional information and engineering support at 1.800.227.1066 or visit us at:

www.danahermotion.com/DMM
www.doverinstrument.com



Helping you build a better machine, faster.